Buy One, Get Two for Free: Simultaneous Knee T2 Mapping and Morphological Analysis On Synthetic Images Using GRAPPATINI

Marion Roux, Tom Hilbert, Jean-Baptise Ledoux, Fabio Becce, Tobias Kober, Patrick Omoumi

Department of Radiology, University Hospital (CHUV), Lausanne,
 Advanced Clinical Imaging Technology (HC CMEA SUI DI BM PI),
 Siemens Healthcare AG, Lausanne
 LTS5, École Polytechnique Fédérale de Lausanne, Lausanne

Unil

UNIL | Université de Lausanne

واللاعدالة مراكر مراكر ومراكل ومراكل مراكل و

## Background

- Most commonly used contrasts to assess morphological joint abnormalities are [1]:
  - Intermediate-weighted TSE (IW)
  - T2-weighted TSE (T2w)
- Quantitative T2 mapping gained interest recently
  - Clinical Relevance: quantitative T2 is a good biomarker of early degenerative disease of cartilage and meniscus [2-4]
- The acquisition of both morphological and quantitative sequences is however time consuming, hindering their spread in clinical use.

Huang M et al. Radiology 2014;273:S1-22.
 Zarins ZA et al. Osteoarthritis Cartilage 2010;18:1408-16.

[2] Kijowski R, et al. Radiology 2013;267:503-13.[4] Rauscher I et al. Radiology 2008;249:591-600.

### **The GRAPPATINI sequence**

### Standard sequences



#### **GRAPPATINI** sequence



## Objective

- Validation of :
  - -(1) T2 values
    - phantom experiment
  - (2) image quality
    - quantitative and qualitative assessment
    - 5 healthy volunteers
    - in comparison to the consecutive acquisition of morphological and quantitative sequences

## 1) Phantom experiments – Methods

- Validation of the T2 estimation with GRAPPATINI
  - Phantom: Tubes filled with different concentrations of agar and gadolinium
  - Experiment: Same hardware and protocols (MESE and GRAPPATINI) used as in in-vivo experiments
  - Reference: Low resolution single slice, single spin-echo sequence



Single Slice Single Spin-Echo (Goldstandard)

T1 Map (MP2RAGE)

Grappatini Knee Protocol

MCSE Knee Protocol

### 1) Phantom Experiments – Results

 MESE > GRAPPATINI > SSSE (stimulated echoes)



### 2) In vivo experiments – Methods

- Single-center prospective study, performed from September 2015 to January 2016
- 5 healthy volunteers (3 males, age 30.2±3.3 years)
- 3T (MAGNETOM Skyra, Siemens Healthcare, Germany), 15channel knee coil
- Acquisition of the **GRAPPATINI** sequence
  - Generation of additional synthetic contrasts with TE=34ms and TE=80ms
- Acquisition of standard morphological TSE images
  - IW (TE=34ms) and T2w (TE=80ms)

## 2) In vivo – Methods – MRI protocol

- 3T (MAGNETOM Skyra, Siemens Healthcare, Germany)
- 15-channel knee coil
- All sequences were acquired in the sagittal plane.

Sequence	Resolution	Number of Slices	TR	TE	Fat Suppression	ТА
GRAPPATINI	0.4x0.4x3mm <sup>3</sup>	36	4880 ms	34 / 80 ms	On	6:22 min
					Total Acq. Time:	6:22 min
IW - TSE	0.4x0.4x3mm <sup>3</sup>	36	3790 ms	34 ms	On	2:52 min
T2w - TSE	0.4x0.4x3mm <sup>3</sup>	36	3790 ms	80 ms	On	2:52 min
MESE	0.5x0.5x3mm <sup>3</sup>	36	1630 ms	∆TE 13 ms	Off	7:09 min
					Total Acq. Time:	12:13 min

# 2) In vivo – Methods – Analysis

- Comparison of synthetic morphological images versus conventional TSE images
  - Quantitative analysis : SNR, CNR
  - Qualitative analysis (2 radiologists blinded to employed sequence, Five-grade scale):
    - Global image quality
    - Image quality of each of the following anatomical structures: cartilage, menisci, cruciate ligaments, bone marrow, muscle, joint fluid, quadricipital and patellar tendons
    - Visual image contrast
    - Visual noise
    - Artifacts

- -2: first image significantly worse than second,
- -1: moderately worse,
- 0: no difference,
- +1: moderately better,
- +2: significantly better



Synthetic Image (TE=34 ms)

Standard IW TSE (TE=34 ms)

## 2) In Vivo Experiments – Results

### Quantitative Analysis

No significant difference in SNR and CNR

### — **SNR** :

 average SNR = 9.9 for both synthetic and conventional sequences, p=0.99

- CNR :

- cartilage/fluid: 6.2 vs. 6.6, p=0.62
- meniscus/fluid: 11.3 vs. 11.6, p=0.81

### 2) In Vivo Experiments – Results

### Qualitative Analysis

- No difference in global image quality or of any of the anatomical structures that were evaluated
  - average score of 0, 95%Cl=[0; 0.4]
- Artefact scores were slightly higher for the synthetic sequences
  - average of -0.1, 95%CI=[-0.002; -0.6]
- Visual noise and contrast were slightly better for the synthetic sequences
  - average score of 0.1 for both, 95%CI=[0.002; 0.6]

### Limitations

- Number of included patients (volunteers)
- Implication in diagnosis imaging

### Conclusion

- T2 map, IW and T2w images can be acquired with one GRAPPATINI sequence in 6.22 min instead of 12.13 min (sum of the acquisition times obtained with the standard techniques)
- GRAPPATINI provides...
  - Accurate T2 values
  - Synthetic IW and T2w images quantitatively and qualitatively similar to conventional TSE images

Buy One, Get Two for Free: Simultaneous Knee T2 Mapping and Morphological Analysis On Synthetic Images Using GRAPPATINI

Marion Roux, Tom Hilbert, Jean-Baptise Ledoux, Fabio Becce, Tobias Kober, Patrick Omoumi

Department of Radiology, University Hospital (CHUV), Lausanne,
 Advanced Clinical Imaging Technology (HC CMEA SUI DI BM PI),
 Siemens Healthcare AG, Lausanne

3 - LTS5, École Polytechnique Fédérale de Lausanne, Lausanne

Canton de Vaudo

UNIL | Université de Lausanne



## **Thank You**

## The GRAPPATINI sequence

### • T2 Map Acquisition:

10-fold undersampled multi-echo spin-echo sequence (MESE)

#### • T2 Map Reconstruction:

GRAPPATINI [7]: Subsequent application of parallel imaging (GRAPPA)
 [5] and model-based reconstruction (MARTINI) [6]

### • Synthetic T2 weighted Images:

 Applying the forward signal Model to simulate TSE contrast with arbitrary echo-time TE

[5] Griswold, Mark A., et al. MRM 47.6 (2002): 1202-1210.
[6] Sumpf, Tilman J., et al. JMRI 34.2 (2011): 420-428.
[7] Hilbert, Tom, et al. ISMRM 22.4077 (2014).